

# Nanodrop\_data\_management

Fay

2023-03-26

## Libraries

```
# install libraries
```

```
library(dplyr)
```

```
## Warning: package 'dplyr' was built under R version 4.2.1
```

```
##
```

```
## Attaching package: 'dplyr'
```

```
## The following objects are masked from 'package:stats':
```

```
##
```

```
##      filter, lag
```

```
## The following objects are masked from 'package:base':
```

```
##
```

```
##      intersect, setdiff, setequal, union
```

```
library(XML)
```

```
## Warning: package 'XML' was built under R version 4.2.3
```

```
library(methods)
```

```
library(plyr)
```

```
## Warning: package 'plyr' was built under R version 4.2.3
```

```
## -----
```

```
## You have loaded plyr after dplyr - this is likely to cause problems.
```

```
## If you need functions from both plyr and dplyr, please load plyr first, then dplyr:
```

```
## library(plyr); library(dplyr)
```

```
## -----
```

```
##
```

```
## Attaching package: 'plyr'
```

```
## The following objects are masked from 'package:dplyr':  
##  
##   arrange, count, desc, failwith, id, mutate, rename, summarise,  
##   summarize
```

```
library(readr)
```

```
## Warning: package 'readr' was built under R version 4.2.1
```

```
library(ggplot2)
```

## Read the tsv nanodrop and write them as csv

### David

```
#read tsv table  
David <- read_tsv("~/GitHub/Namibia_project/Data/Nanodrop_measurements/David_Nanodrop_DNA-Extractions.tsv")  
  
# remove ffirst column  
David <- David[, -1]  
  
#add column sample type  
David <- David %>%  
  mutate(animal = "David")  
  
#change the column names  
  
write.csv(David,  
  "~/GitHub/Namibia_project/Data/Nanodrop_measurements/CSV/David.csv",  
  row.names = FALSE)
```

### Düppel

```
#read tsv table  
Duppel_1 <- read_tsv("~/GitHub/Namibia_project/Data/Nanodrop_measurements/Düppel_26022023_2.tsv")  
  
# remove ffirst column  
Duppel_1 <- Duppel_1[, -1]  
  
#change the column names  
  
write.csv(Duppel_1,  
  "~/GitHub/Namibia_project/Data/Nanodrop_measurements/CSV/Düppel_26022023.csv",  
  row.names = FALSE)
```

### Rodents

```

#read tsv table
Rodent_1 <- read_tsv("~/GitHub/Namibia_project/Data/Nanodrop_measurements/Rodents_26032023.tsv")
Rodent_2 <- read_tsv("~/GitHub/Namibia_project/Data/Nanodrop_measurements/Rodents_27032023.tsv")

#combine the measuerements
Rodent <- rbind(Rodent_1,Rodent_2)

# remove first column
Rodent <- Rodent[,-1]

#change the column names

write.csv(Rodent_1,
          "~/GitHub/Namibia_project/Data/Nanodrop_measurements/CSV/Rodents_26032023.csv",
          row.names = FALSE)

```

## Merging the tables

```

Nanodrop <- rbind(Duppel_1, Rodent)

# remove spaces
colnames(Nanodrop) <- gsub(" ", "_", colnames(Nanodrop))

Nanodrop <- Nanodrop %>%
  dplyr::rename(Quality_260_280 = '260/280',
                Quality_260_230 = '260/230')

```


## Assessing the quality of the DNA

```

ggplot(Nanodrop, aes(x = Nucleic_Acid, y = Quality_260_280)) +
  geom_jitter() +
  labs(x = "Concentration of DNA in ng/µl", y = "Quality of DNA 260/280 ratio",
       title = "Assessing quality of DNA with Nanodrop")

```


## Assessing quality of DNA with Nanodrop



Let's filter out the samples with extremely high ratios

```
ggplot(Nanodrop %>%
  filter(Quality_260_280 < 2.5, Quality_260_280 > 1.5 ),
  aes(x = Nucleic_Acid, y = Quality_260_280)) +
  geom_jitter() +
  labs(x = "Concentration of DNA in ng/μl", y = "Quality of DNA",
    title = "Assessing quality of DNA with Nanodrop")
```


## Assessing quality of DNA with Nanodrop



Assesing the quality of 260/230


```
ggplot(Nanodrop,  
  aes(x = Nucleic_Acid, y = Quality_260_230)) +  
  geom_jitter() +  
  labs(x = "Concentration of DNA in ng/μl", y = "Quality of DNA, 260/230",  
    title = "Assessing quality of DNA with Nanodrop")
```

### Assessing quality of DNA with Nanodrop




```
ggplot(Nanodrop %>% filter(Quality_260_230 < 60),  
  aes(x = Nucleic_Acid, y = Quality_260_230)) +  
  geom_jitter() +  
  labs(x = "Concentration of DNA in ng/μl", y = "Quality of DNA, 260/230",  
    title = "Assessing quality of DNA with Nanodrop")
```

## Assessing quality of DNA with Nanodrop



```
ggplot(Nanodrop %>%  
  filter(Quality_260_230 < 2.5, Quality_260_230 > 1.8),  
  aes(x = Nucleic_Acid, y = Quality_260_230)) +  
  geom_jitter() +  
  labs(x = "Concentration of DNA in ng/μl", y = "Quality of DNA, 260/230",  
       title = "Assessing quality of DNA with Nanodrop")
```

## Assessing quality of DNA with Nanodrop



Select the samples having the golden ratios

```
Nanodrop %>%
  filter( Quality_260_280 < 2.5, Quality_260_280 > 1.5 )
```

```
## # A tibble: 104 x 11
##   Sampl~1 User_~2 Date_~3 Nucle~4 Unit A260_~5 A280_~6 Quali~7 Quali~8 Sampl~9
##   <chr>   <chr>   <chr>   <dbl> <chr>   <dbl>   <dbl>   <dbl>   <dbl>   <chr>
## 1 S01    fay_w    26/03/~  1     "ng/~  0.019  0.008   2.38   0.33 DNA
## 2 S11    fay_w    26/03/~  3.1   "ng/~  0.062  0.038   1.63   0.24 DNA
## 3 S12    fay_w    26/03/~  10.3  "ng/~  0.206  0.108   1.92   0.34 DNA
## 4 S28    fay_w    26/03/~  1.8   "ng/~  0.037  0.017   2.16   0.03 DNA
## 5 S42    fay_w    26/03/~  30.9  "ng/~  0.618  0.336   1.84   1.84 DNA
## 6 BW24   fay_w    26/03/~  18.7  "ng/~  0.374  0.201   1.86   0.78 DNA
## 7 BW40   fay_w    26/03/~  63.4  "ng/~  1.27   0.68    1.87   1.13 DNA
## 8 BW44   fay_w    26/03/~  11.3  "ng/~  0.226  0.112   2.01   0.21 DNA
## 9 C01    fay_w    26/03/~  4.4   "ng/~  0.088  0.044   2      -1.03 DNA
## 10 C03   fay_w    26/03/~  6.8   "ng/~  0.135  0.088   1.54   1.81 DNA
## # ... with 94 more rows, 1 more variable: Factor <dbl>, and abbreviated
## #   variable names 1: Sample_ID, 2: User_name, 3: Date_and_Time,
## #   4: Nucleic_Acid, 5: `A260_(Abs)`, 6: `A280_(Abs)`, 7: Quality_260_280,
## #   8: Quality_260_230, 9: Sample_Type
```